

Gloves and Glove Dispensers



Acrylic Glove Dispensers

Clear, durable acrylic dispenser keeps workspaces neat. Wipes clean. Shipped so that they can be mounted on a wall if necessary.

Size

Quad, 21" H x 4-1/4" W x 11-1/2" L
(53.34 cm x 10.8 cm x 29.21)

Double, 10" H x 4-1/4" W x 11-1/2" L
(25.4 cm x 10.8 cm x 29.21)

Triple, 16" H x 4-1/4" W x 11-1/4" L
(40.64 cm x 10.8 cm x 28.5)

Single Horizontal, 6" H x 4-1/2" W x 11-1/2" L
(15.24 cm x 11.43 cm x 29.21)

Triple Narrow, 30" H x 5-1/2" W x 4" L
(76.2 cm x 13.97 cm x 10.16)

Single Vertical, 10" H x 5-1/2" W x 4-1/4" L
(25.4 cm x 13.97 cm x 10.8)

Double Narrow, 20" H x 5-1/2" W x 4-1/4" L
(50.8 cm x 13.97 cm x 10.8)

Product No.

BA070

BA071

BA072

BA073

BA074

BA075

BA076



Stainless Steel Glove

Individual interlocking stainless steel rings provide a remarkable sense of touch while protecting the hand from cuts. Designed to be worn between two latex gloves, the stainless steel glove provides protection from infection. The Velcro® wristband holds the glove comfortably in place. Universal glove - fits left or right hand. UL Listed.

Size

X-small

Small

Medium

Large

X-large

Product No.

BE050

BE052

BE054

BE056

BE058

Have questions? **Our Customer Service Department is friendly AND knowledgeable**