

Scales



Easy-to-Use Digital Weight Indicator

LCD backlit readout makes for excellent display in all types of lighting. Indicator and scale can be easily washed down. This durable unit can be wall, bench or tower mounted. Indicator has a RS232 bi-directional interface hookup for use with optional printer.

Stabilization Time: 1-4 seconds
 Weighing Modes: kilograms, grams, ounces, pounds
 Power Requirements: AC Adapter (included), 120V or 6 alkaline "C" batteries or NiCad
 Platform 12" x 12" x 4" (30.48 cm x 30.48 cm x 10.16 cm)
 Constructed of stainless steel

All Stainless Steel Washdown Scale Base

Rugged High-Resolution Scale for the Surgical Pathology Laboratory & Morgue.

Totally constructed of 100% stainless steel including all nuts, bolts and washers. Durable and completely washable. Includes AC Adapter.

Capacity: 50 lbs. (22.7 kg)
 Readability: .005 LB/.002 kg
 Platform Size: 12" x 12" (30.5 cm x 30.5 cm)

Product No. BB027

Portable Postmortem Scale Stand

The portable stand brings the scale to a convenient location and helps reduce reaching and lifting of postmortem organs. The offset center post reduces the possibility of tipping over.

The compact and stable platform stand is completely fabricated of stainless steel for easy cleaning and durability. Four easy rolling, swivel and locking casters provide complete mobility. For use with BB002 and BB014 scales.



Base is Stainless Steel



Heavy-Duty Locking Casters

*BB035 Scale Stand
(Scale Sold Separately)*

Dimensions: 24" L x 24" W x 78.8" H
 (60.96 cm x 60.96 cm x 200.15 cm)

Product No. BB035