

ACCULAB Scales



BB020

Pocket Pro® PP401

This balance is the ultimate in portability and the perfect choice for sample-prep. Featuring a back-lit display. Simple operation and a unique pan design which protects the scale during transportation and storage.

Product No. BB020



BB021

Electronic EC-2100

Designed to handle rough laboratory use. Balance weighs in grams or ounces. The handy hold button allows reading to remain on display after removing specimen. Large plastic platform and angled LCD are easy to use and read. Scale is battery operated and offers full calibration and accuracy over a wide weight range.

Product No. BB021



BB022

VIC-4kg

Large sturdy stainless steel platform, bright angled LCD, keypad calibration, tare and memory functions. AC or battery operation (AC adapter included).

Product No. BB022

Acculab Scales	PP401 Product No. BB020	EC-2100 Product No. BB021	VIC-4kg Product No. BB022
Capacity:	400 g	2100 g	4100 g
Readability:	0.1 g	1 g	1 g
Linearity:	+/-2d	±1 g	+ (-) 1d
Weighting Modes:	grams, ounces	grams, ounces, pounds	grams, ounces, pounds
Calibration:	Auto Calibration	Auto Calibration	Auto Calibration
Tare:	To Capacity	To Capacity	To Capacity
Auto Off:	Yes	Yes	Yes
Platform Size:	2.76" x 2.95" (7.0 cm x 7.4 cm) Plastic	5.3" x 5.5" (13.4 cm x 13.9 cm)	5.5" x 5" (13.97 cm x 12.7 cm) Stainless Steel
Power:	AA	9V	9V Alkaline, 9V Nicad AC Adapter (included)*
Dimensions:	2.95" x 4.33" x .8" (7.4 cm x 10.9 cm x 2.0 cm)	8.5" x 6" x 2" 21.5 cm x 15.2 cm x 5.0 cm	6.8" x 3.15" x 9.8" (17.2 cm x 8.0 cm x 24.8 cm)
Weight:	1 lb.	1.5 lbs..	2.9 lbs.
Warranty:	12 Months	12 Months	12 Months
Optional Bowl (Not shown BB024): Heavy-Duty Round Plastic Bowl allows for easy clean-up. Holds 1050 ml			