

## ACCULAB Scales



BB010



BB060

### Portable Scale

The VIC-Series portable top-loader features auto calibration, tare, four weighing modes, audio function tone, large stainless platform with a built-in lock down bracket, large LCD display, spill rings and easy to use "up front" controls. Unit operates on AC adapter (included) or 9V battery. One year warranty included.

For VIC-10 kg Dimensions: 13 cm sq. order BB010

For VIC-6 kg Dimensions: 13 cm sq. order BB060

Portable Scale	VIC-10 kg	VIC-6 kg
	<b>Product No. BB010</b>	<b>Product No. BB060</b>
Capacity:	10100 g/352.7 oz	6100 g/211.6 oz.
Weighting Modes:	grams and ounces	grams and ounces
Calibration:	Auto Calibration from Keypad	Auto Calibration from Keypad
Tare Range:	To Capacity by Subtraction	To Capacity by Subtraction
Overload Indication:	OL on LCD	OL on LCD
Underload Indication:	UL on LCD	UL on LCD
Auto Off:	2, 5 and 10 min.	2, 5 and 10 min.
Readability	1 g	1 g
Linearity	+ (-) 2d	+ (-) 2d
Dimensions:	6.8" x 3.15" x 9.8" (17.2 cm x 8.0 cm x 24.8 cm)	
Platform Size:	5.5" x 5" (13.97 cm x 12.7 cm)	