

Digital Scales



Compact Electronic Balance Scale

Simple and easy to operate, this rugged laboratory scale is designed with a high impact plastic housing and a large detachable plastic (or optional stainless steel) weighing pan.

With full digital calibration and accuracy over a wide temperature range, this scale is battery operated. Optional AC Adapter available.

Capacity:	2000 g
Display:	LCD .31" H (8 mm)
Readability:	1 g
Linearity:	+/- 2 g
Weighing Modes:	grams, oz, lb, kg
Tare:	To Capacity
Calibration:	Full Digital Calibration with Front Key
Auto Off:	5 min. (Separate Re-zero Function)
Pan Size:	5.35" W x 5.35" D (136 mm x 136 mm)
Power:	6 AA Batteries (not included) or AC Adapter (optional)
Dimensions:	(5.51" W x 7.79" D x 1.75" H) (140 mm x 198 mm x 44.5 mm)
Weight:	Approx. 1.1 lb. (455 g)
Warranty:	12 Months

Product No. BB026

Optional Adapter Product No. BB036



Compact Digital Scale

This rugged, splash proof digital scale is designed for rough laboratory use. Balance weighs in grams or ounces. Large 1" LCD display for easy accurate reading. Large removable stainless steel or plastic weighing pan. Battery operated for portable use. Optional AC adapter available.

Capacity:	2000 g
Display:	7 Segment LCD, 1" H (25 mm)
Readability:	1.0 g or .0021 lbs.
Linearity:	+/- 2 g
Weighing Modes:	grams, oz, kg, lb
Tare:	To Capacity
Calibration:	Factory Set (Re-Zero Function)
Auto Off:	5 min. Battery Only
Pan Size:	9.05" W x 7.48" D (230 mm x 190 mm)
Power:	6 D Batteries (not included) or AC Adapter (optional) Low Battery Indicator
Dimensions:	9.61" W x 9.13" D x 5.39" H (244 mm x 232 mm x 137 mm)
Weight:	Approx. 3.53 lbs. (1.6 kg)
Warranty:	12 Months

Product No. BB033

Optional Adapter Product No. BB036