


# Laboratory Dissecting Lights

## Portable Dissecting Light

The BB108 model has a single lighthead mounted on a floor stand with casters. The 26" (66 cm) extension arm rotates 60° around the mounting post. Yoke rotates 180° around the arm and the lighthead 180° within the yoke. The top of the post is 5'-1" (155 cm) high. The lighthead can be raised to 6' (183 cm) and lowered to 3'-1" (94 cm).

**Product No. BB108**



*Portable Dissecting Light  
BB108*



## Single Track Ceiling Mount Light

Mopec's most recommended dissection light. Vertical adjustment of 48" (122 cm) with 360° lighthead rotation. Available with one lighthead on a single track and with two lightheads (not shown) one on each of the two tracks. Carriage travel within the track is 48" (122 cm). Vertical centerline of the lighthead is approximately 38" (96.5 cm) from the vertical centerline of the suspension tube when the two arm sections are fully extended horizontally. This provides a nominal 26" (66 cm) extension beyond the side arm edge of the typical autopsy table. Optional intensity control is available for both models.

**Single Lighthead Track Light - Product No. BB112**

**Dual Lighthead Track Light - Product No. BB113**