

Laboratory Dissecting Lights

Dissecting Lighting II™

- 2000 Footcandles at 24" (61 cm)
- Color temperature of 3600° K
- Three 50 watt quartz halogen bulbs Wired in parallel to assure continuous illumination; 2000 hour bulb life
- Blue filters and dichroic coatings keep beam cool
- Prefocused spot
- Spot size of 13" (33 cm) diameter
- Drift-free, self-balancing arm; shadow-free design
- Option for drop tube lengths for ceiling mount - choice of 30" or 42"
- Includes SteriHandle™ and disposable light cover sample
- 115V Standard
- UL Listed, CSA certified
- 1-Year warranty



Ceiling or Wall-Mounted Dissecting Light with Single Lighthead

Lightweight – inexpensive lighting. Single lighthead mounted on the ceiling suspension or wall bracket. These features are the same as BB113. The ceiling mount has 360° continuous rotation around the suspension tube. The wall mount moves parallel with the surface it's mounted on.



Mounting Bracket - BB111

Ceiling-Mounted Light -Product No. BB110

Wall-Mounted Light - Product No. BB111



Ceiling-Mounted Dissecting Light with Dual Lightheads

The BB109 model has two lightheads mounted on a single ceiling suspension. The lightheads can be rotated 250° about the suspension tube and has a 19" (48 cm) vertical travel range upward from horizontal.

Further rotation can be achieved through the second arm and yoke. The lighthead itself rotates 180° on its own axis. The vertical center of the lighthead is approximately 53" (134 cm) from the suspension tube when the two arms and yoke are fully extended. Maximum lighting using minimum space.

Product No. BB109