

X-Ray Illuminators




X-Ray Illuminator
One Viewing Area
Desk or Surface Mount - BB115
Recessed Mounting - BB118

Deluxe X-Ray Illuminators


- Baked-on white enamel finish
- Thin compact design
- 3-Lamp light source. Viewers operate on 118 V 60 Hz and are supplied with a three conductor, line cord. Special voltages are available.
- Low leakage trigger start ballast
- Surface and recessed models
- UL listed and CSA approved

The deluxe illuminators are only 3-3/8" (8.6 cm) deep and are ideal for installations that require minimum projection. The deluxe illuminators feature one piece construction for each tier with steel channel reinforcement of the rear section. The viewing area is lighted by three fluorescent tubes and is controlled by its own individual on/off rocker switch. Optional automatic film activated switches are also available.

Each illuminator is supplied with a gravity grip file holder which extends the entire width of the viewer. A one-piece plexiglass viewing panel is also supplied for each tier. Illuminators can be ordered as a desk or surface mount (portable) or as a recessed mount (permanent).



X-Ray Illuminator
Two Viewing Areas
Desk or Surface Mount - BB116
Recessed Mounting - BB119



X-Ray Illuminator
Three Viewing Areas
Desk or Surface Mount - BB117
Recessed Mounting - BB120



Completed lab

Desk or Surface Mount, One Viewing Area
14" W x 21" H x 3-3/8" D
(35.6 cm x 53.3 cm x 8.6 cm)

BB115

Recessed Mounting, One Viewing Area
14-3/4" W x 21-3/4" H x 3-3/8" D
(37.5 cm x 55.3 cm x 8.6 cm)

BB118

Desk or Surface Mount, Two Viewing Areas
28" W x 21" H x 3-3/8" D
(71.1 cm x 53.3 cm x 8.6 cm)

BB116

Recessed Mounting, Two Viewing Areas
28-3/4" W x 21-3/4" H x 3-3/8" D
(73 cm x 55.3 cm x 8.6 cm)

BB119

Desk or Surface Mount, Three Viewing Areas
42" W x 21" H x 3-3/8" D
(106.7 cm x 53.3 cm x 8.6 cm) (Not Shown)

BB117

Recessed Mounting, Three Viewing Areas
42-3/4" W x 21-3/4" H x 3-3/8" D
(108.6 cm x 55.3 cm x 8.6 cm)

BB120