

Portable Supply Cart and Podium



Heavy-Duty Stainless Steel Construction



Front view

All autopsy supplies, instruments forms and other presentation materials can be stored on two sturdy 24" x 24" fixed shelves. Four 5" locking casters allow the unit to be easily moved throughout the autopsy suite. The large writing surface has a lip for containment of any size forms, or requisition and the writing surface height starts at 43" and angles up to a height of 53". Overall dimensions: 53" tall x 24" wide x 24" deep.

Custom shapes and styles are available.

Product No. BJ100



Adjustable Height Portable Writing Podium

Provide a clean, convenient and organized place to collect and analyze your data. The sloped design provides a spacious 25" long x 20" deep writing surface. A clear acrylic cover on the back panel will protect and secure the document for viewing. The podium has a 19" wide x 19" deep locking metal drawer for storing supplies. The unit can be adjusted in height to provide additional comfort. The four locking casters provide easy mobility and will lock in place when needed.

Product No. BJ101

