

Staining Dishes



Coplin Glass Staining Jar with Glass Cover

Heavy-duty molded construction with glass cover. Ten slide grooves accommodate 20 - 3" x 1" slides (back to back). Inside Dimensions: 3.35 L 1.75" W x 2.44" W x 1.38" D (85 mm x 62 mm x 35 mm).

Product No. BA001

Coplin Glass Staining Jar with Glass Cover

Heavy, clear glass jar with clear glass cover. Accommodates 5 3" x 1" slides (back to back). Dimensions: 4.75 H 1.75" W x 3.75" D (120.65 mm x 44.45 mm x 95.25 mm).

Product No. BA002

Coplin Glass Staining Jar with Plastic Cover

Heavy, clear glass jar, 5 molded internal grooves, tight fitting screw-on-plastic black cover prevents fluid evaporation. Accommodates 10 3" X 1" slides (back-to-back). Dimensions: 4.3" H x 1.6" Dia. (109 mm x 40.64 mm).

Product No. BA003

Vertical Wide-Top Heavy Clear Glass Jar

Jar has 8 slots, each slot holding 2 3" x 1" slides (back-to-back). The wide top is convenient for transferring slides from jar to jar. Glass cover included. Dimensions: 4.375" x 2.25" x 1.5". (111 mm x 57 mm x 38 mm).

Product No. BA004