

Slide Storage



BG006
60 Slide
Staining Dish

BG002
30 Slide
Staining Dish

BG008
60 Slide
Staining Tray

BG009
5 Slide
Staining Tray

BG004
30 Slide
Staining Tray

Stainless Steel Staining Dishes

Stainless steel construction with covers

30 Slide Capacity, Dimensions: 4.5" L x 3.5" W x 3" H

Product No. BG002

60 Slide Capacity, Dimensions: 8" L x 3.75" W x 3" D

Product No. BG006

Glass Slide Staining Dish

Heavy-duty molded construction with glass cover. 10 slide grooves accommodates 20 - 3" x 1" slides (back-to-back). Inside dimensions: 3.35" L x 2.44" W x 1.38" D (85 mm x 62 mm x 35 mm).

Product No. BG001

Stainless Steel Staining Trays

Removable handle, stainless steel construction, accepts 3" W x 1", 1-1/2" or 2" H slides.

30 Slide Capacity

Product No. BG004

60 Slide Capacity

Product No. BG008

10 Slide Capacity

Product No. BG009

(attached handle)