

Slide Storage



6 Drawer Slide Storage Cabinet shown with base cabinet

Microscope Slide Storage Cabinet

Stackable cabinets are made from heavy-duty sheet metal with a baked-on enamel finish. Can be locked together for stability. Drawers are removable for easy handling. Knobs are rounded for quick pull. Each drawer comes with a tab for labeling. Holds 5,000 slides!

6 Drawer Slide Storage Cabinet

15-3/4" W x 18-3/4" D x 5" H
(40 cm x 47.62 cm x 127 cm)

Product No. BH011

Cabinet Drawers for Slide Cabinet

15-3/4" W x 18-3/4" x 2-3/4"
(40 cm x 47.62 cm x 6.9 cm)

Product No. BH012



BH014

Reinforced Cardboard
Slide Tray

BH008

Plastic Slide Folder

Reinforced Cardboard Slide Tray

Two doors provide protection in a dust-free environment for safety while carrying or storing slides. Holds 20 standard 3" x 1" slides. Each slide is raised off the surface allowing easy removal simply by pressing on the end. Slides are placed in individual compartments. Closed Dimensions: 7-1/2" W x 11-1/2" L (19.05 cm x 29.21 cm).

Product No. BH014

Plastic Slide Folder

Made of durable white plastic, folder holds 20 standard 3" x 1" slides. Transparent cover. Each slide rests on four plastic points to prevent sticking, and provides easy removal. Dimensions: 7-1/2" W x 1/4" D x 11-1/2" H (19.05 cm x .63 cm x 29.21 cm).

Product No. BH008