

Replacement Fume Exhaust Filters



Activated Charcoal Bulk

Activated Charcoal Bulk

For use with tissue processors. 3-1/4" lb. container.

Order No. BF004



Charcoal Filter - BF002

Charcoal Filter

Manufactured with nearly five pounds of coconut shell-derived, disposable activated charcoal material encased in a heavy-gauge perforated metal frame.

For use with Mopec's BF Series Fume Hoods.


12" W x 12" L x 1" H

Order No. BF002

For use with Mopec's BF400 and BF500 Grossing Workstations.

7" W x 14" L x 1" H (not shown)

Order No. BF008



Particle Filter

Particle Filter

Constructed of high efficiency polyester fiber encased in a heavy cardboard frame, disposable filters will trap and retain most particulate matter larger than 5 microns. For use with Mopec's BF Series Fume Hoods.

11-7/8" W x 11-7/8" L x 1/2" H

Order No. BF003



Potassium Permanganate Filter

Potassium Permanganate Filters

Safely absorbs and neutralizes formaldehyde fumes with powerful potassium permanganate disposable filters.

Use with Mopec's MB and ME Series Grossing Workstations. For models prior to 2004. Dim: 6-1/4" W x 14" L x 1" H

Order No. BF001

For use with Mopec's BF Series Fume Hoods.

Dim: 12" W x 12" L x 1" H (not shown)

Order No. BF005

For use with Mopec's BF400 and BF500 Grossing Workstations.

Dim: 7" W x 14" L x 1" H

Order No. BF007

Use with Mopec's newer MB and ME Series Grossing Workstations.

Dim: 4" W x 15-3/4" L x 3" H

Order No. BF009



Potassium Permanganate Filter for MB & ME Models