

MB100 Grossing Station / Counter Top

Counter top Grossing Workstation with Small Specimen Sink



MB100 shown with optional MB006 Leg Frame

Keep The Workflow Moving

Ideal for small laboratories with limited space constraints, the MB100 compact countertop unit may also serve as a secondary station in larger, busier facilities requiring additional workspace to meet increased demands. The MB100 workstation will maximize work space with standard conveniences such as recessed fluorescent lighting, a small specimen sink with wrist blade handles and swing spout faucet.

The versatile MB100 can also be converted into a free standing unit with the optional stainless steel leg frame. Our optional heavy-duty frame features four leveling feet for secure placement on uneven surfaces.

STANDARD FEATURES:

- A** Washable In/Cm Ruler
- B** Polyethylene Cutting Board (Not Shown)
- C** G.F.C.I. Outlet
- D** Magnetic Instrument Tool Bar
- E** Fluorescent Light
- F** S.S. Paper Towel Holder
- G** Perforated Vent Grill
- H** Faucet w/Swing Spout
- I** Small Specimen Sink

OPTIONAL FEATURES:

See MB Options page

