

# JD500 Cadaver Tray Lift Bars

## Cadaver Lift Bar with Straps



The JD500 Cadaver Lift Bar System product consists of two heavy-duty bars along with nylon straps that securely attach to the hand slots of the body trays for easy elevation. To be used in conjunction with the JD Series Cadaver Lifts. Made of stainless steel with a #4 satin finish. All welds are a continuous type and smooth, free of burrs and brushed to remove any discolorations.

### JD500 Standard Features


**Support Bars:** 1-1/2", 11 gauge wall

**Nylon Straps:** 1", 36 inch basket


**Basket Style:** Rated at 2,200 lbs.

**Fork Saddle:** 14 gauge stainless steel

**Product No. JD500**



*JD500 shown with JE100 Storage Mattress*



*JD500 shown with JC100 Body Tray*