

JD600 Cadaver Lift Body Straps

Cadaver Lift Body Bar with Straps



*JD600 Cadaver Lift Bar
shown with JD710 Body Lift*

The Mopec cadaver lift body straps are designed in conjunction with the JD-Series Cadaver Lifts. Simply slide the nylon straps under the body and attach to the top elevating bar for easy and effortless cadaver removal. Made of stainless steel with a #4 satin finish. All welds are a continuous type and are smooth, free of burrs and brushed to remove any discoloration. Lifting capacity is 1,000 pounds.

JD600 Standard Features

Support Bars: 1-1/2 inch, 11 gauge wall

Nylon Straps: 3 inch, 54 inch basket

Fork Saddle: 14 gauge stainless steel

Product No. JD600

Check out

MOPEC CARRIES A WIDE ASSORTMENT OF BODY BAGS, INCLUDING THE HEAVY-DUTY BODY BAG WITH SIDE PANEL.

See Pages
95-97



JD600 Cadaver Lift Bar shown with DD100 Autopsy Cart