

# JD1000 Smart Body Lift – A Mopec Exclusive


"No Hand Cranking Needed"



JD1000 shown with JC100 Body Tray and JD400 Roller Pallet. Activation button on flexible cord allows lifting from various positions

mobility in confined spaces and narrow aisles, the maximum capacity of the JD1000 is 750 lbs. The elimination of constant cranking and low maintenance make this the perfect lift for any lab or medical center. Delivered fully assembled.

The JD1000 is a compact, battery powered cart that helps alleviate physical stress caused by manual transporting and lifting. Easy to operate, the JD1000 has a built-in charger. Charging is accomplished with a single 110 volt outlet. Depending on the desired height and load, a single charge will easily lift hundreds of times. Battery recharge time (total discharge) is 4 hours. Lightweight and highly maneuverable with easy



JD1000 shown with JC100 Body Tray

## TECHNICAL DATA

- **Max lifting weight capacity:** 750 lbs.
- **Overall height:** 83"
- **Base width:** 36"
- **Nest height (up position):** 78"

## ERGO FORCE DATA

### Two hand push from:

- **Rear (empty):** 6.0 - 7.5 ft/lbs.
- **Rear (300 lb . load):** 10.5 -12.75 ft/lbs
- **Side (empty):** 4.0 - 6.5 ft/lbs
- **Side (300 lb load):** 5/5 - 6.75 ft/lbs

## LIFT SPEED AND TRAVEL

- **Up/Down Travel:** 52"
- **Speed (empty):** 3.5" per second
- **Speed (300 lb. load):** 2.7" per second

**Product No. JD1000**