

## JE100 Body Tray Mattress

Stainless Steel



Cadaver Lift shown with JD500 Lift Bars with Straps and JE100 Body Mattress

### Safe and Easy Cadaver Storage

Made of stainless steel, this unique mattress was designed to fit all standard body trays and transport carts. Once placed under the cadaver, four sturdy eyebolts on each corner of the mattress connect easily to any body lift or transport system.

The mattress also allows elevation of the cadaver during autopsy for fluid drainage and cleaning of the work surface. Mopec Body Lift sold separately.

### Storage Mattress

Description	Product No.
Body Tray Mattress	JE100



JE100 Storage Mattress shown with JD500 Lift Bars with Straps and JC100 Body Tray



JD710 Cadaver Lift Shown with JD500 Lift Bar with Straps and JE100 Storage Mattress

**Check out**

[www.mopec](http://www.mopec.com)

FOR ONLINE SPECIALS