

CD400 Fold-Up Wall Mount Table



*CD400 Fully opened.
One hand easy operation*

When you don't have room for a complete workstation, install the CD400, which provides a compact work area for multiple procedures. Prep, surgery and exams can be performed on this table. Designed to be wall mounted, this table drops to a full 90°. When finished, simply fold and lock in an upright position against the wall.

Made of 18 gauge stainless steel, the rim has a marine edge around perimeter which will contain fluids on the work surface. Load capacity is 300 lbs. Table dimensions are 20" W x 45" L.

CD400 Standard Features

- 18 Gauge stainless steel with all exposed welds ground and polished.
- Axle assembly uses 2 nylon wheels that allow the table to be stored in an upright, vertical position when not in use.
- Stainless steel release handle unlocks the table from an upright position to a 90° horizontal position.
- Four holes are positioned on table for mounting



CD400 Folds up or down to 90°.



CD400 Folded