

Dissection Aids



Binocular Magnifying Loupe

German Optics 1.8X magnification with coated optics. Lightweight, mounted on spectacle frame. 30 cm focal distance, 10 cm viewing area, pupillary distance adjustment. Complete with carrying case.

Product No. AK100



Professional Padded Instrument Case

Now you can keep your instruments neat and protected. Zippered Professional Instrument Case is made of black padded vinyl and lined with a red brushed velvet-like cloth. Elasticized safety straps (38) are attached to securely hold instruments. Protective safety flap is attached to side. Measures 25" (635 mm) from end to end when fully opened. 12" W x 16" H (305mm x 406.4 mm) when zippered shut. Outside of case can be wiped down with a damp cloth. **Instruments not included.**

Product No. AA310

Student Instrument Case

Designed for the student to carry instruments back and forth. Zippered Student Instrument Case is made of vinyl and lined with a dark red brushed velvet-like cloth. Twelve elasticized safety straps are attached to securely hold instruments. Pen strap is in the middle of case. Protective safety flap is attached to side. Measures 15" (381 mm) from end to end when fully opened. 7" W x 8-1/4" H (178 mm x 209.5 mm) when zippered shut. Outside of case can be wiped down with a damp cloth. Choose either burgundy or black. **Instruments not included.**

Product No. AA300 - Black case (Not Shown)

Product No. AA301 - Burgundy case

