

Rulers



Precision Caliper


Excellent for measuring all size tissue specimens. Made of stainless steel. Sliding caliper with a finger tip nut to lock down ruler. Total length is 4" (101.6 mm) with graduated etched markings in inches and millimeters. Width capacity is 3.5" (88.9 mm).

Product No. BA040

8" Stainless Steel Foldable Ruler

Tired of using your ruler and coming up short? Our stainless steel ruler is a full 8 feet long. It also folds in half and has a 1/4 diameter hole on either end for easy storage and can be hung on the wall. Total length is 8' (243.84 cm) with plastic faced ruled edge in inches and centimeters. Width capacity is 1" (2.54 cm).

Product No. BA048



Did you Know **we carry PathAceto Wipes?**

Disposable acetic acid cloth mordant wipes

See page 119

Product No. BE037