

# Bone Dust Collector

The HEPA Bone Dust Collector has all the same features as our QuietCut vacuum. Designed to attach to Mopec, Stryker and many other autopsy saws, the vacuum motor is wired to start when the saw is turned on.

## Powerful Vacuum System:

Harmful airborne pathogens are trapped in a uniquely designed disposable HEPA filter.

## Mobile

The Bone Dust Collector is mounted onto a single lightweight module that can easily roll from the autopsy suite to the surgical dissection bench.



**HEPA Filter Replacement**  
Product No. BD053

*BD080 For use with Stryker Saw (220V Saw Not Available)*

## HEPA Vac Specifications U.S. Specifications

Electrical	115V/60Hz - 8.5A
Cord Length	10 ft. (305 cm)
Dimensions	23" L x 15" W x 23-1/2" H (58.4 CM X 38.1 X 59.7cm)
Weight	25 lbs (11.3 kg)
Hose	1" I.D. x 102" (2.5 cm x 259 cm)

## Bone Dust Collector

Electrical:	115V/50-60Hz, 8.5 Amps
Electrical Cord Length:	120"
Hose Length:	1.25" I.D. x 120" (3.2 cm x 30.5 cm)

When ordering, be sure to inform Mopec's Customer Service what style autopsy saw will be used.

HEPA Bone Dust Collector with Mopec Saw Adapter  
**Product No. BD075**

Mopec 220V - For Export **Product No. BD076**

HEPA Bone Dust Collector with Stryker  
**Saw Adapter Product No. BD080**

Stryker 220V - For Export **Product No. BD081**

