

Storage Containers, Strainers, Funnels



Standard Bowl Strainer

Used for tissue or specimen storage. Made of stainless steel.

Size

1-1/4 qt. (1.2 liter) (4-1/8" Body Dia.,
5-3/4" D x 4-3/4" Over Bead Dia.)
(10.5 cm x 14.6 cm x 12.1 cm)

2 qt. (1.9 liter) (4-7/8" Body Dia.,
6-3/4" D x 5-1/2" Over Bead Dia.)
(12.4 cm x 17.2 cm x 14.0 cm)

4-1/4 qt. (4 liter) (6-1/8" Body Dia.,
7-5/8" D x 7-1/8" Over Bead Dia.)
(15.6 cm x 19.4 cm x 18.1 cm)

Cover, Fits BA018

Cover, Fits BA019

Cover, Fits BA020

Product No.

BA018

BA019

BA020

BA021

BA022

BA023



Standard Bowl Strainer

Mopec's standard stainless steel bowl strainer offers a broad collection bin. Small scattered drainage holes (app. 300 holes - 1/4" in diameter) allow slow controlled straining and rinsing. Capacity 14 quart. Excellent for small delicate tissue specimens.

Size

16-3/16" Dia. x 5-1/2" Deep
41.11cm x 13.97 cm.)

Product No.

BA050



Funnels

Size - Can hold

2-3/8 oz. (68 ml)

13 oz. (369 ml)

2 qt. (1.9L)

Dia. of Funnel

3"

5"

8-3/4"

Product No.

BA042

BA043

BA044