

Utility Bin



Multi-Purpose Specimen/ Instrument Utility Bin with Strainer

Constructed of rugged washable polyethylene, this perforated multi-purpose bin with strainer provides an ideal unit for temporary specimen fixation, instrument soaking and treatment or specimen collection and presentation. Approximately 300 holes, 1/4" diameter. Side handles allow unit to easily lift from the base for fast efficient drainage and/or rinsing. Unit includes a depressed molded lid to minimize exposure to potentially harmful fixation or disinfectant solutions. Stackable.

Size:

Top, 20-1/2" L x 16" W (52.07 cm x 40.64 x 12.70 cm)

Strainer, 20" L x 15" W x 5" D (50.8 cm x 38.10 cm x 12.70 cm)

Base, 20" L x 15-1/2" W x 7" D (50.8 cm x 39.37 cm x 17.78 cm)

Product No. BA053

Mopec can supply you
with all your Autopsy Accessories

Dissecting Specimen Pins

Securely holds specimen with pins. Large quantity package of 100 pins. Size: 1-1/4" (3.81 cm).

Product No. AD010



Use for specimens - not an actual specimen



Actual size

