


Polyethylene Dissecting Boards



Mopec's unique photo blue dissecting board provides an excellent photography background



Dissecting Boards with Legs

Elevate your dissecting board

Mopec offers legs on both our polyethylene and neoprene cork dissecting boards. Legs screw into board with a variety of leg lengths available.

Size	Product No.	
	White	Blue
11" W x 12" L x 3/4" Thick w/ 3" Legs (28 cm x 30.5 cm x 2 cm)	BC003	n/a
23" W x 16" L x 3/4" Thick (58.4 cm x 40.6 cm x 1.9 cm) BC001 with 3" (7.6 cm) Legs	BC021	BC034
16" W x 12" L x 3/4" Thick (40.6 cm x 30.5 cm x 1.9 cm) BC002 with 3" (7.6 cm) Legs	BC022	BC035

Polyethylene Dissecting Boards

Manufactured of durable 3/4" (1.9 cm) thick, easy-to-clean polyethylene with rounded edges and non-skid rubber bumpers (non-skid leg tips shown above on Model BC003). Opaque white surface provides excellent background for cutting surgical or autopsy tissue. Mopec's unique photo blue dissecting board provides an excellent photography background.

Size	Product No.	
	White	Blue
23" W x 16" L x 3/4" Thick (58.4 cm x 40.6 cm x 1.9 cm)	BC001	BC031
16" W x 12" L x 3/4" Thick (40.6 cm x 30.5 cm x 1.9 cm)	BC002	BC032
Ruler for Dissecting Boards (Cut to fit, with adhesive back)	BC029	

**Oversized...
Undersized...
Most sizes available...
Custom Sizes are No Problem!**

**Mopec can
assist
you**

