



# Dissecting Tables



BC006 Fixed Leg Table  
(Shown with Dissecting Board BC031)



BC010 Drain Tray

12" (304.8 mm) Wide lip



BC010 Drain Tray  
(Shown with Dissecting Board BC031)

## Stainless Steel Dissecting Table with Fixed Legs

The table top is tapered toward the middle drain for fast and easy drainage. A drain hose can be connected to a fitting on the underside of the tray for control of fluid. The sturdy 12" (30.48 cm) legs can be set on any counter or autopsy table and accommodates large or small dissection boards.

Tray Dimension: 25" L x 18" W x 1" D  
(63.5 cm x 45.7 cm x 2.5 cm)

Overall Height: 12" (30.5 cm)  
Stainless Steel Table with Fixed Legs

**BC006**

## Stainless Steel Dissecting Drain Tray

A multi-purpose unit, this stainless steel tray can be used to confine spills or contamination to a restricted area. All of Mopec's dissecting boards will fit into the one inch deep tray. The sloped end overhangs any sink, allowing all fluids to drain directly into sink.

28" L x 18" W x 1" D  
(71.12 cm x 45.7 cm x 2.54 cm)

Dissecting Drain Tray

**BC010**



*Elevate your dissecting board with Mopec's portable, all stainless steel dissecting tables*