

Body Bags



BE001
Lightweight Center
Zipper Style

BE002
Heavy-Duty
with Handles

BE003
Heavy-Duty
with Handles
Envelope
Zipper Style

BE004
Lightweight
Envelope
Zipper Style

Lightweight Center Zipper Style

Material	8 ml. Virgin Plastic
Size	36" x 93" (91.4 cm x 236.2 cm)
Girth	60" (152.4 cm)
Thread	Nylon
Zipper	Heavy Duty #10
Weight	1 lb. 10 oz. (.7 kg)
Capacity	Over 200 lbs. (91 kg)
Handles	No

BE001

Heavy-Duty with Handles Center Zipper Style

Material	Polyester
Size	36" x 93" (91.4 cm x 236.2 cm)
Girth	60" (152.4 cm)
Thread	#69 Nylon
Zipper	#5 Nylon, 88" (223.5 cm)
Weight	6 lbs. (.8 kg)
Capacity	At least 300 lbs. (136 kg)
Handles	Six 2" (5.1 cm) Nylon Webbing Runs Under Bag

BE002

Heavy-Duty with Handles Envelope Zipper Style

Material	Polyethylene
Size	34" x 93" (86.36 cm x 236.2 cm)
Girth	60" (152.4 cm)
Thread	#69 Nylon
Zipper	#5 Nylon, 140" (355 cm) with Two Tabs Can be Opened at Any Point
Weight	6 lbs. (2.7 kg)
Capacity	At least 300 lbs. (136 kg)
Handles	Six 2" (5.1 cm) Nylon Webbing Runs Under Bag

BE003

Lightweight Envelope Zipper Style

Material	8 ml. Virgin Plastic
Size	36" x 93" (91.4 cm x 236.22 cm)
Girth	60" (152.4 cm)
Thread	Nylon
Zipper	Heavy Duty #10
Weight	1 lb. 10 oz. (.7 kg)
Capacity	Over 200 lbs. (91 kg)
Handles	No

BE004