

Bodyrests



Autopsy Head Drape

Minimize infection and contamination exposure during autopsy by containing particulate matter.

- Easily assembled and disassembled
- Adjusts easily over torso or head
- Just cut side slots to allow for easy access during autopsy

Autopsy Head Drape (Frame, includes 4 head-drapes)

Product No. BE010

Autopsy Head Drape Only (10 per Package)

Product No. BE011



Slide Body Transfer Board

Perimeter handles offer safe alternative to moving cadavers from cooler to morgue table

Ideal for transferring cadavers between multi-level work surfaces

Provides secure handling before and after autopsy procedures

Wash with any standard hospital disinfectant.

Dimensions: 72" L x 17" W x 1/4" D
(182.88 cm x 43.18 cm x .63 cm)

Product No. FD010



Front View



Support/Elevating Body Slats

Set of four body slats elevate the body on a table during the autopsy or morgue preparation. While elevated, waste drains without obstruction. Suggested placement of slats are shoulders, buttocks, heels and head. Made of a chemically resistant white plastic, they each measure 3" wide (7.62 cm) and 22-5/8" long (57.5 cm). Can be cleaned using water and one of Mopec's specially formulated cleaners.

Set of Four

Product No. BC011