



# Pathology Dissecting Knife Holder

*Safely Organize Your Knives*



This simple design knife holder allows the handle to protrude for safe gripping. Simply place your sharp knives in either slot and they are kept out of harm's way.

The two tubes make it easy to organize your instruments and the angled notches keep the knives separated, ready for use. The slot at the bottom allows for liquid drainage and easy cleaning of the holder.

Made of a heavy-duty, durable plastic ensures a long product life and reduces dulling of your dissecting knives.




**Mopec The Better Product People**  
**800-362-8491 or [www.mopec.com](http://www.mopec.com)**

## WHAT SIZE KNIFE DO YOU HAVE?

- ✓ 14 " (35.56cm) Large Section Knife
- ✓ 12 " (30.48cm) Brain Sectioning Knife
- ✓ 10 " (25.40cm) Lung Knife
- ✓ 8 " (20.32cm) Autopsy Knife
- ✓ 4 " (10.16cm) Cartilage Knife
- ✓ Scalpel blade handles

Almost any size or shape knife used in your pathology lab can be stored

Product Number AH116



*For a complete catalog detailing our products and services, please call or visit our web site*



21750 Coolidge Highway • Oak Park, MI 48237  
Toll Free: 800-362-8491 • Tel: 248-291-2040 • Fax: 248-291-2050

[www.mopec.com](http://www.mopec.com)