

## OE100 Dualdraft Workstation



### Countertop Dualdraft Workstation

Ideal for small laboratories with limited space constraints that require a high demand for air ventilation. The OE100 provides a compact, trimming environment featuring our unique dual-draft ventilation system. As the worksurface accumulates trimming peripherals and tissue, the downdraft and backdraft will automatically adjust to the changes. Airflow is constant thus air velocity will increase accordingly in the open area.

Mopec offers an exceptionally wide variety of standard features. Provided with conveniences such as recessed fluorescent lighting, hot/cold water faucet, wrist blade handles and adjustable shelving. The versatile OE100 is designed to be placed on an existing counter top, or free standing with the optional stainless steel leg frame. The heavy-duty frame features four leveling feet for secure placement on uneven surfaces.

## OE100 Dualdraft Workstation

### ME100 STANDARD FEATURES:

- A Fixture w/Swing Spout
- B Paper Towel Holder
- C Adjustable Shelf
- D Fluorescent Light (Under Canopy)
- E G.F.C.I. Outlet
- F Perforated Grid Plates
- G Backdraft Grill
- H Downdraft Chamber

### ME100 OPTIONAL FEATURES:

- 00001 Self Contained Exhaust System - Single
- 00003 Formalin Dispensing Container
- 00004 Dictation Equipment Stand
- 00006 Leg Frame
- 00007 Foot Pedal
- 00008 Perimeter Rinse
- 00009 Safety Splash Shield
- 00010 Hand Spray Rinse
- 00011 Magnifier for Halogen Light
- 00014 Video Camera Arm
- 00015 Halogen Task Light Dual
- 00016 Camera Stand
- 00018 Formalin Collection System
- 00019 Literature Dispenser
- 00032 Computer CPU Bracket
- 00033 LCD Flat Screen Monitor/Keyboard Adjustable Arm System
- 00035 Shelves and Adjustable Brackets

