

# Specimen Transport and Storage Systems

## Pathport® 3 Mobile Transport and Specimen Storage

New for 2007

# PATHPORT<sup>3</sup>

MOBILE SPECIMEN STORAGE

Pathport<sup>3</sup> is a stainless steel constructed transport and storage system, specifically manufactured to accommodate pathology specimens. Equipped with a custom fit, disposable formaldehyde-neutralizing pad, the Pathport<sup>3</sup> is specially designed to absorb harmful spills and vapors. Spill clean-up can be accomplished by simply replacing the pad.

For full O.S.H.A. compliance, order the Pathport Transport System, exclusively from Mopec.

**Dimensions:** 43"H x 18"W x 36"L  
(108cm H x 45cm W x 92cm L)

Formalin and other pathology fixatives are classified as regulated waste. O.S.H.A.'s Bloodborne Standard 29 CFR Part 1910.1030 mandates that pathology specimens be placed in secondary containers which are:

- **Constructed to contain all contents and prevent leakage of fluids during handling, storage, transport or shipping**
- **Labeled**
- **Closes prior to removal to prevent spillage of contents during handling, storage, transport and shipping.**

## Specimen Transport and Storage All in One Product!



*A better alternative for pathology specimen transport and storage. Pathport<sup>3</sup> offers a safe, secure, and economical method for transporting and storing specimens.*

**Convenience:** Offers the combined convenience of secure specimen storage and transportation anywhere.

**Safe:** Pathport<sup>3</sup> allows the transport and storage of pathology specimens in full compliance with O.S.H.A. bloodborne and hazard communication standards.

**Secure:** Pathport<sup>3</sup> offers a secure alternative to store any pathology specimen.

**Effective:** Pathport<sup>3</sup> disposable pads effectively absorb and neutralize harmful formalin spills within the Pathport<sup>3</sup> Transport System.

**Economical:** Pathport eliminates costly annual exposure monitoring and hazard communication training for transport personnel and laboratory staff.

**Product No. BE020**

**Disposable Replacement Pads**

(1 per package)

36" L x 18" W x 1/4" Thick  
(90 cm x 45 cm x .7 cm)

**Product No. BE021**