

MB Series Optional Features

Customized to Meet Any Laboratory Specifications

safely collect residual formalin for safe disposal or recycling. Formalin not included.

MB020 Literature/Form Holder

Lexan literature dispenser capable of holding 1" letter-size documents for easy handling and referral (three tier).

MB026 Ventilated Trash Receptacle

Removable ventilated trash container, 12 gallon capacity. Included with a clear flexible ventilation duct, which is tied into the grossing station's ventilation system.

MB028 Formalin Dispensing System

Dispensing system that is contained within the grossing station. Includes 5-1/2 gallon dispensing carboy, pump and backsplash mounted valve. Indicator lights for high/low formalin liquid levels in carboy. Formalin not included.

MB029 Computer CPU Bracket

Adjustable size mounting bracket for CPU, keyed to mounting plate with two pins for easy removal or installation. Computer not included.

MB031 Removable Writing Ledge

Add an additional writing ledge surface to your grossing station with the twist of a knob. Just slide it over the sink and you have an extended writing or working ledge that is easily removable.

MB041 Combination Ledge and Pencil Drawer

Convenient pull out writing ledge and utility drawer made of stainless steel construction.

MB047 Pneumatic Foot Pedal

Pneumatic foot operated water control. Foot operated toggle on/off switch. Hand operated mixing valve next to faucet.

MB048 Stainless Steel Shelf

Stainless steel shelf adjusts to accommodate storage needs.

MB049 New LCD Flat Screen Monitor & Keyboard Adjustable Arm System

Allows positioning of LCD display for greater visibility Up, down, forward and back. Full range of motion for maximum comfort level keyboard adjustable arm. Effortlessly position a keyboard in the most comfortable typing position, seated or standing use. Flat screen, monitor, keyboard and computer not included.

MB050 Ventilated Trash Receptacle on Wheels

Removable ventilated trash container, 12 gallon capacity with 3" swivel casters with brake mechanism. Included with a clear flexible ventilation duct, which is tied into the grossing station's ventilation system.

BL800 Waste Disposal

1 HP commercial grade disposal is air-cooled for continuous operation. Power switch activates solenoid to continuously feed water to the unit. Stainless steel internal parts prevent corrosion. Swivel impellers ensure the disposal will never jam.

Mopec has developed a complete line of grossing workstations to fit your individual laboratory requirements. Each workstation can be customized with a variety of available options—all at an affordable price.



MB015

MB016

MB017

MB019

MB020

MB026



MB031

MB041

MB047

MB048

MB049

MB050

BL800