

## MB Series Optional Features

### Customized to Meet Any Laboratory Specifications

safely collect residual formalin for safe disposal or recycling. Formalin not included.

#### **MB020 Literature/Form Holder**

Lexan literature dispenser capable of holding 1" letter-size documents for easy handling and referral (three tier).

#### **MB026 Ventilated Trash Receptacle**

Removable ventilated trash container, 12 gallon capacity. Included with a clear flexible ventilation duct, which is tied into the grossing station's ventilation system.

#### **MB028 Formalin Dispensing System**

Dispensing system that is contained within the grossing station. Includes 5-1/2 gallon dispensing carboy, pump and backsplash mounted valve. Indicator lights for high/low formalin liquid levels in carboy. Formalin not included.

#### **MB029 Computer CPU Bracket**

Adjustable size mounting bracket for CPU, keyed to mounting plate with two pins for easy removal or installation. Computer not included.

#### **MB031 Removable Writing Ledge**

Add an additional writing ledge surface to your grossing station with the twist of a knob. Just slide it over the sink and you have an extended writing or working ledge that is easily removable.

#### **MB041 Combination Ledge and Pencil Drawer**

Convenient pull out writing ledge and utility drawer made of stainless steel construction.

#### **MB047 Pneumatic Foot Pedal**

Pneumatic foot operated water control. Foot operated toggle on/off switch. Hand operated mixing valve next to faucet.

#### **MB048 Stainless Steel Shelf**

Stainless steel shelf adjusts to accommodate storage needs.

#### **MB049 New LCD Flat Screen Monitor & Keyboard Adjustable Arm System**

Allows positioning of LCD display for greater visibility Up, down, forward and back. Full range of motion for maximum comfort level keyboard adjustable arm. Effortlessly position a keyboard in the most comfortable typing position, seated or standing use. Flat screen, monitor, keyboard and computer not included.

#### **MB050 Ventilated Trash Receptacle on Wheels**

Removable ventilated trash container, 12 gallon capacity with 3" swivel casters with brake mechanism. Included with a clear flexible ventilation duct, which is tied into the grossing station's ventilation system.

#### **BL800 Waste Disposal**

1 HP commercial grade disposal is air-cooled for continuous operation. Power switch activates solenoid to continuously feed water to the unit. Stainless steel internal parts prevent corrosion. Swivel impellers ensure the disposal will never jam.

*Mopec has developed a complete line of grossing workstations to fit your individual laboratory requirements. Each workstation can be customized with a variety of available options—all at an affordable price.*



**MB015**

**MB016**

**MB017**

**MB019**

**MB020**

**MB026**



**MB031**

**MB041**

**MB047**

**MB048**

**MB049**

**MB050**

**BL800**