

# CC100 Autopsy Sink

Autopsy Sink Wall Mount - Center Cart Approach 114"



CO016 Optional Scale Stand Mounted on CC100

CC100 Autopsy Sink shown with DA100 Autopsy Cart

The ergonomic design of our CC100 Autopsy Sink eliminates the physical strain of moving bodies. The autopsy can be performed directly on the cart.



Optional CO020 Organ Rinse Basket



Autopsy Cart Latch



Dissection Area Grid Plates


## CC100 Autopsy Sink

### Autopsy Sink Wall Mount - Center 114"

The CC100 center approach wall mounted autopsy sink provides all of the features normally found in a typical pedestal type autopsy table without the space requirements. The autopsy cart overhangs the autopsy station allowing fluids to drain directly into the sink. On each side of the sink are dissection areas of equal size with removable and interchangeable grid plates.

The dissection area has grid plates which are elevated above a sloped surface that drains directly into the sink. Wall mounted design allows for easy floor cleaning.

**Product No. CC100**



NOTE: Sink sprays are shown mounted on skirt panel

#### STANDARD FEATURES:

- A** Perforated Grid Plates (4)
- B** Autopsy Cart Latch
- C** Fixture w/Swing Spout (2)
- D** Body Rinse Fixture (2)
- E** Large Instrument Drawer (2)
- F** Sink Spray Assy. (2)
- G** GFCI Outlet (2)
- H** Sink Well
- Rinse Basket (Not Shown)

#### OPTIONAL FEATURES:

- I** BL300 1/2 HP Disposal
- CO002 Autopsy Sink Leg Frame (Not Shown)
- CO003 Perimeter Rinse
- M** CO006 Eye Wash
- K** CO007 Exhaust Duct and **N** Ventilation Grill (2)
- CO008 Fluorescent Light Assembly (2)
- CO016 Scale Stand Wall Mounted (Not Shown)
- CO020 Organ Rinse Basket (Not Shown)
- CO022 Evidence Screen (Not Shown)
- K** CO006 Eye Wash

