

Autopsy Saw and Saw Reconditioning



Stryker Autopsy Saw

The Stryker autopsy saw has long been a reliable instrument for the pathology industry. The high speed oscillating action quickly cuts through bone without harming soft tissue. Used for removing the cranial cap, making linear cuts or sectioning small bone specimens. Durable and easy to handle. UL listed.

Product No. BD001



Saw Reconditioning Service

Mopec can recondition your old autopsy saw, often at a fraction of the price of a new saw.

Mopec will replace:

- All bearings
- All worn internal parts
- Electrical switch and power cord
- Armature and electrical fields

An evaluation fee will be charged if the saw cannot be reconditioned.

6-Month Unconditional Warranty!

NOTE: This applies to Mopec saws only. For Stryker Autopsy Saws, please contact Stryker. Call our toll free number today, (1-800-362-8491)

Product No. BD009

Did you know **you can order online 24/7?**