

# Neoprene Cork and Disposable Boards



## Neoprene Cork Dissecting Boards

These dissecting boards consist of a 1/4" layer of neoprene cork permanently cemented to a 3/4" marine plywood base. Neoprene cork provides a very durable cutting surface resistant to absorption while permitting thorough rinsing.

### Size

23" W x 16" L x 1" Thick  
(58.4 cm x 40.6 cm x 2.5 cm)

### Product No.

**BC004**

16" W x 12" L x 1" Thick  
(40.6 cm x 30.5 cm x 2.5 cm)

**BC005**

## Neoprene Cork Boards with Legs

See Legs On Next Page.

### Size

BC004 with 3" (7.6 cm) legs

### Product No.

**BC023**

BC005 with 3" (7.6 cm) legs

**BC024**



## Disposable Cutting Boards

When contamination control is important

- Foam inner liner for pinning tissue
- Standard and metric scale on cutting surface
- Reference areas for assessing specimen size
- Ideal for specimen photography
- Excellent for teaching institutions
- Ideal for infectious tissue cases

### Size

Small, 48 per Package  
6-1/4" x 7-1/2" (15.87 cm x 20.9 cm)

### Product No.

**BC007**

Medium, 24 per Package  
6-1/4" x 15-3/8" (15.87 cm x 39.05 cm)

**BC008**

Large, 12 per Package  
13-1/8" x 15-3/8" (33.33 cm x 39.05 cm)

**BC009**