

# MB200 Grossing Station / Counter Top



The MB200 countertop model is ideal for smaller labs or used as a secondary station for larger facilities. Also available as a free standing unit



BL300  
Heavy-Duty Disposal




MB019  
Formalin Collection



MB004  
Dictation Equip. Stand



MB011  
Safety Splash Shield



MB017  
Digital Camera Stand

## MB200 Grossing Station / Counter Top

### Counter Top Grossing Workstation with Large Sink




### The MB200 Features One of the Deepest Sinks in the Industry

Designed for small laboratories with limited space constraints, the MB200 compact countertop unit may also serve as a secondary station in larger, busier facilities requiring additional workspace to meet increased demands. The MB200 workstation will maximize work space with standard conveniences such as recessed fluorescent lighting, wrist blade handles and a swing spout faucet.

Also incorporated is a large integrally constructed sink (19" x 14" x 10" deep), allowing for easy and convenient handling of large specimens and/or containers.

The versatile MB200 may also be converted into a free standing unit with the optional stainless steel leg frame. Our optional heavy duty frame features four leveling feet for secure placement on uneven surfaces.

MB200 shown with optional MB006 Leg Frame



#### STANDARD FEATURES:

- A** Washable In/Cm Ruler
- B** Polyethylene Cutting Board (Not Shown)
- C** G.F.C.I. Outlet
- D** Magnetic Instrument Tool Bar
- E** Fluorescent Light (Under Canopy)
- F** S.S. Paper Towel Holder
- G** Faucet w/Swing Spout
- H** Sink Well
- I** Perforated Vent Grill

#### OPTIONAL FEATURES:

See MB options page