

# MB400 Grossing Station




Mopec's MB400 free standing workstation features more usable space and an extra deep sink



MB039/MB040  
Stainless Steel Organizer  
Bin Dispensers



MB016 Dual Halogen  
Lights with MB013  
Magnifying Lens



MB017 Digital Camera  
Stand



MB008 Flammable Storage Cabinet. Cabinet  
features two doors, an adjustable shelf and a key  
lock for added security

## MB400 Grossing Station

### Free Standing Floor Model with Large Sink



Designed for any size laboratory, the MB400 includes all the conveniences and space to perform the most demanding grossing procedures. Twelve inches longer than the MB200, the MB400 incorporates ample work space along with standard conveniences such as a recessed fluorescent light, large specimen sink, effective ventilation and ergonomic location of all amenities.

The MB400 also incorporates the large integrally constructed sink (19" x 14" x 10" deep), allowing for easy handling of large specimens and/or containers. Included with a base cabinet and frame, the MB400 grossing station is free standing and can be placed virtually anywhere in the lab. The envelope of the unit is sized according to standard architectural parameters and will fit into most cabinet schemes.

**NOTE:** If MB008 - Flammable storage option is chosen, you cannot add drawers. Option MB008 adds an extra inch to entire height of grossing station.

#### STANDARD FEATURES:

- A** Washable In/Cm Ruler
- B** Polyethylene Cutting Board (Not Shown)
- C** G.F.C.I. Outlet
- D** Adjustable Shelf
- E** Fan Control Module (Only when MB005 option is selected)
- F** Magnetic Instrument Tool Bar
- G** Fluorescent Light
- H** S.S. Paper Towel Holder
- I** Perforated Vent Grill
- J** Support Frame
- K** Support Cabinet
- L** Faucet w/Swing Spout
- M** Sink Well
- N** Stainless Steel Organizer Bins

#### OPTIONAL FEATURES:

See MB options page

