

# MB600 Elevating Grossing Station


One of Our Newest and Most Advanced Grossing Workstations



MB039/MB040 Stainless Steel Organizer Bin Dispensers



MB005 Exhaust System



MB026 Trash Ventilated Receptacle



MB009 Perimeter Rinse System


## MB600 Elevating Grossing Station

### Elevating Grossing Workstation with Extra Large Sink



The MB600 is the most technically advanced backdraft grossing station available. Incorporating all of the standard features of the MB100, MB200 and MB400, the adjustable height MB600 workstation meets all the requirements of any laboratory.

The MB600 includes an open workspace allowing ample room for those who work standing or sitting. This elevating station creates an ergonomically correct dissection area for any physician or technician. The MB600 is also wheelchair accessible.



The MB600 features Mopec's Backdraft ventilation

#### STANDARD FEATURES:

- A** Magnetic Instrument Tool Bar
- B** Polyethylene Cutting Board
- C** Washable In/Cm Ruler
- D** Fan Control Module (Only when MB005 Option Is Selected)
- E** Paper Towel Holder
- F** Hot/Cold Water Fixture
- G** Elevating Up/Down Switch
- H** G.F.C.I. Duplex Outlet
- I** Fluorescent Light
- J** Sink Well
- K** Stainless Steel Organizer Bin (Small and Large)

#### OPTIONAL FEATURES:

See MB options page

