

MB700 Elevating Wall Mounted Grossing Station

Our Most Advanced Grossing Workstation Just Got Better!



The MB700 allows for easy clean-up without obstruction



MB005 Exhaust System



MB019 Formalin Collection



MB009 Perimeter Rinse System


MB700 Elevating Wall Mounted Grossing Station

Our Most Advanced Grossing Workstation Just Got Better!

A worldwide leader in innovative equipment, the MB700 is a Mopec exclusive. Save valuable space with the MB700, the most technically advanced backdraft wall mounted grossing station available. The MB700 elevates from 29 inches to 41 inches. Incorporating all of the standard features of the MB600, the MB700 meets all the requirements of any laboratory. Additional options can be added making this model a multi-functional unit. Floors can be thoroughly cleaned without obstruction. The MB700 includes an open workspace allowing ample room for those who work standing or sitting. This elevating station creates an ergonomically correct dissection area for all laboratory personnel. The MB700 is also wheelchair accessible.



Elevates from 29" to 41" to accommodate all technicians



STANDARD FEATURES:

- A** Magnetic Instrument Tool Bar
- B** Washable In/Cm Ruler
- C** Paper Towel Holder
- D** Hot/Cold Water Fixture
- E** Elevating Up/Down Switch
- F** GFCI Duplex Outlet
- G** Fluorescent Light
- H** Sink Well
- I** Stainless Steel Organizer Bins (Small and Large)

OPTIONAL FEATURES:

See MB options page