

HH100 Handicapped Workstation



HH100 shown with optional features.

HH100 Handicapped Accessible Grossing Dissection Workstation

HH100 handicapped grossing workstation has all the convenient features needed for gross dissection and the large open leg space gives ample room to accommodate a wheelchair. Stainless steel sink is standard and is 16" x 12" x 4" deep. This multi-functional dissection table is designed with all fixtures within easy reach for the technician.

Standard Features

Full Length Storage Shelf: 10" deep x 64" long

Hot/Cold Water Fixture: Swing spout faucet with built in vacuum breaker, has easily accessible wrist blade handles for hot and cold water control.

Electrical Receptacle: Receptacle is GFCI type with watertight cover.

Optional Features

Disposal: BL800 1 HP commercial grade disposal


Hand Spray Rinse: Soft spray head for easy positioning.

Unit can be customized to fit your lab requirements with other additional options. Contact Customer Service with your specific needs.

HH100 Handicapped Workstation



If you have special requirements for your grossing area, use the Mopec Solution. We will be happy to design and manufacture the perfect workstation for you



STANDARD FEATURES:

- A** Hot/Cold Water Fixture
- B** Stainless Steel Shelf
- C** GFCI Outlet
- D** Vent Tie-in
- E** Sink Well

OPTIONAL FEATURES:

- F** B1800 - 1 HP Disposal
- G** Hand Spray Rinse