

# HI100 Dissection Table – Vented Right Hand Sink HI150 Dissection Table – Vented Left Hand Sink



## Grossing/Dissection Table with Ventilation Hood

**Excellent for Autopsy Suites:** The HI100 is exclusively designed to connect to any facility ventilation system. This versatile and economical dissection table with ventilation hood offers an efficient and maintenance free ventilated workstation. The ample work surface is created toward the large sink for fast drainage. The optional disposal can be utilized for complete sanitary removal of residual tissue. The dissection table incorporates many standard design features.

### Standard Features

**Work Surface:** Standard design feature includes a recessed working area which has a built-in positive slope for drainage. The slope is 3/4" and is pitched

toward the large sink located on the left side of the unit. Table top with double flange construction is designed for heavy-duty usage. Material utilized is 14 gauge stainless steel with a #4 satin finish.

**Water Fixture:** Hot/Cold water fixture with single lever mixing valve is conveniently located on the center of the large sink section. The spout includes a vacuum breaker for water protection.

**Electrical:** A GFCI protected duplex receptacle is included, completely wired and ready for use.

**Frame:** The frame is constructed of 1-1/2" square tubing, 16 gauge wall. All welded construction with each tube having a continuous bead, weld free of burrs and discoloration. Adjustable stainless steel feet included for leveling.

**Vent Hood:** Fabricated of 18 gauge stainless steel with a #4 satin finish. The ventilation grill is 8" high at

## HI100 Dissection Table – Vented Right Hand Sink HI150 Dissection Table – Vented Left Hand Sink

an angle of 60 degrees. The ventilation grill is perforated with 3/16" x 5/8" slots on 13/16" x 3/8" staggered centers, yielding a 30% open area.

**Construction:** The superior construction utilizes TIG welding on all seams and mating sections which are then ground smooth and regrained to match the #4 satin factory finish.

**Plumbing:** All internal plumbing and electrical lines are factory installed and require only trade hookup.

### Optional Features

**Disposal:** BL800, Heavy-duty commercial disposal 1 HP, 115/1/60. Complete with on/off switch, water solenoid and vacuum breaker plumbed and wired ready for use.

**Hand Spray Rinse:** Soft spray head for easy positioning provides an additional shield between you and specimen.

Unit can be customized to fit your lab requirements with other additional options. Contact Customer Service with your specific needs.

#### STANDARD FEATURES:

- A** Sink Well
- B** Dissection Area
- C** Ventilation Hood
- D** Vent Tie-in
- E** Water Fixture

#### OPTIONAL FEATURES:

- F** BL800 – 1 HP Disposal
- G** HO007 – Hand Spray Rinse

