

# OA700 Necropsy Table



OA700 modified extra long and with extra grid plates

The OA700 Series Necropsy Tables offer a long ventilated work surface, which is an effective means of pulling the odors away from the user. The built-in baffle system allows consistent airflow from one end of the table to the other. Full sized integrally constructed table available with table flush system, the table can remain clean and sanitary.

## Standard Features

**Base & Plenum:** Fabricated of 18 and 16 gauge stainless steel with a #4 satin finish.

**Hot/Cold Fixture:** Hot and cold water fixture with convenient wrist blade handles; gooseneck faucet and vacuum breaker fixture mounted directly on faucet.

**Perforated Grid Plates:** Three grid plates are fabricated of 14 gauge stainless steel and are completely interchangeable.

**Inner Baffling:** The inner baffling is full length and completely removable for cleaning. Fabricated of 18 and 16 gauge stainless steel with a #4 satin finish.

**Electrical Receptacle:** Two (2 each) electrical receptacles are provided and are located on each side of the pedestal. The receptacles are GFCI type with watertight cover.

**Hand Spray Rinse:** Soft spray head for easy and controlled rinsing of work surface.

## OA700 Necropsy Table

### Optional Features

**BL800:** Heavy-duty 1 HP commercial disposal, complete with on/off switch, water solenoid, and vacuum breaker plumbed and wired ready for use.

**Rinse Bar:** Full rinse virtually eliminating the accumulation of debris and fluids. Control valve included for water flow adjustment.

**Sink Grid Plate:** Constructed of stainless steel, positioned under faucet on low end of unit.

**Scale Stand Table Mounted:** Portable stand support



Optional Sink Grid Plate

#### STANDARD FEATURES:

- A** Three Grid Plates
- B** GFCI Outlet
- C** Hot/Cold Fixture
- D** Damper with Adjustable Locking Knob
- E** Hand Spray Rinse

#### OPTIONAL FEATURES:

- BL800 1 HP Disposal (Not Shown)
- OO008 Rinse Bar (Not Shown)
- OO026 Scale Stand Table Mounted (Not Shown)
- OO029 Sink Grid Plate (Not Shown)

