


Double Edge Blade Handle and Blades

AJ111



AJ127

AJ119

Blade Handle & Replacement Blades  
AJ125



AJ124

AJ112

No. 8 Handle  
For #60/#70 Blades  
AJ115



AJ114

AJ113

## Scalpel Handles and Blades

### Double Edge Blade Handle

Mopec's double edge blade handle is ideal for trimming surgical specimens and is a versatile tool for any laboratory. The double edge blade is securely held in place, yet easy to replace. The pencil style handle provides a firm yet mobile grip. Double edge blades are sold 5 per package.



### Stainless Steel Blade Remover

Made for the safe removal of scalpel blades from handles. Easy to use - just insert handle into the blade remover and while holding the remover's edges down, pull out the handle. It's that quick and easy. Made of stainless steel, the blade remover is also autoclavable.

Product No. AJ107



1. Insert blade and handle into blade remover



2. Holding blade remover firmly down, pull out handle.

A Blade Handle and Blades	
<b>Description</b>	<b>Product No.</b>
Double Edge Blade Holder	AJ127
Double Edge Razor Blades	AJ111

### Blade Handles

B Blade Handles	
<b>Description</b>	<b>Product No.</b>
Blade Handle, Metal Handle	AJ125
Replacement Blades for disposable knife	AJ124
No. 8 Handle for #60/#70 Blades	AJ115
Stainless Steel Handle, (Smooth) for #60/#70 Blades	AJ113
Stainless Steel Handle, (Fluted) for #60/#70 Blades	AJ114
Stainless Steel Handle, (Smooth) for #60/#70 Blades	AJ113

### Dissecting Surgical Blades

C Carbon Steel	
<b>Size</b>	<b>Non-Sterile</b>
Size #60	AJ112
Size #70	AJ119



### Scalpels - Sterile

AJ130

One per box unless otherwise noted.

Scalpels	
Scalpel, 1-1/2" (38.1 mm) Blade	AJ130
Round Tip Blade, Pkg. of 10	AJ124
Blade Handle (metal)	AJ125

### Disposable Surgical Blade Remover

Safely removes dull contaminated blade from surgical handle and reduces the number of accidental cuts. The safe way to remove and dispose of blades when disposal facilities are not always available. Simply insert and press handle to dispense old blade safely into sealed container. Can hold approximately 100 blades. Works with any brand of scalpel handle. Not compatible with double sided blades. 10 blade removers per package.

Product No. AJ116



1. Insert handle and blade in remover



2. Press and hold while pulling out handle



3. Blade is now in container