

OE300 Dualdraft Trimming Workstation

Designed for larger capacity laboratories, the OE300 provides a compact, well ventilated trimming environment featuring our unique dualdraft ventilation system. As the worksurface accumulates trimming peripherals and tissue, the downdraft and backdraft will automatically adjust to changes in open area. Airflow is constant thus air velocity will increase accordingly. The OE300 is designed as a free standing unit and includes a heavy-duty frame featuring four leveling feet for secure placement on uneven surfaces.



Shown with optional shelves

OE300 STANDARD FEATURES:

- A** Fixture w/Swing Spout
- B** Paper Towel Holder
- C** Adjustable Shelf
- D** Fluorescent Light (Under Canopy)
- E** G.F.C.I. Outlet
- F** Perforated Grid Plates
- G** Backdraft Grill
- H** Downdraft Chamber
- I** Leg Frame

OPTIONAL FEATURES:

See me options page.

